

Product Datasheet

Human Chemokine C-C-Motif Receptor 3 (CCR3) ELISA Kit (orb1208131)

Description	Human Chemokine C-C-Motif Receptor 3 (CCR3) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51677
Sensitivity	11.2pg/mL
Entrez	1232
Expiration Date	Please enquire.