

Product Datasheet

Human Clara Cell Protein 16 (CC16) ELISA Kit (orb1208139)

Description	Human Clara Cell Protein 16 (CC16) ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid , cell culture supernates or other biological fluids.
Uniprot ID	P11684
Sensitivity	23.5pg/mL
Entrez	7356
Expiration Date	Please enquire.