

Product Datasheet

Human B-Cell Activation Factor Receptor (BAFFR) ELISA Kit (orb1208163)

Description	Human B-Cell Activation Factor Receptor (BAFFR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96RJ3
Sensitivity	0.061ng/mL
Entrez	115650
Expiration Date	Please enquire.