

Product Datasheet

Human Angiopoietin Like Protein 1 (ANGPTL1) ELISA Kit (orb1208178)

Description	Human Angiopoietin Like Protein 1 (ANGPTL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O95841
Sensitivity	0.055ng/mL
Entrez	9068
Expiration Date	Please enquire.