

Product Datasheet

Human Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit (orb1208186)

Description	Human Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6UWV6
Sensitivity	1.39ng/mL
Entrez	339221
Expiration Date	Please enquire.