

Product Datasheet

Human Alcohol Dehydrogenase 1 (ADH1) ELISA Kit (orb1208194)

Description	Human Alcohol Dehydrogenase 1 (ADH1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P07327
Sensitivity	0.72ng/mL
Entrez	124
Expiration Date	Please enquire.