

Product Datasheet

Human Alpha-1-Microglobulin/Bikunin Precursor (a1M) ELISA Kit (orb1208205)

Description	Human Alpha-1-Microglobulin/Bikunin Precursor (a1M) ELISA Kit
Reactivity	Human
Range	4.68-300ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, urine or other biological fluids.
Uniprot ID	P02760
Sensitivity	1.8ng/mL
Entrez	259
Expiration Date	Please enquire.