

Product Datasheet

Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) ELISA Kit (orb1208212)

Description	Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q4QQW2
Sensitivity	5.2pg/mL
Entrez	171369
Expiration Date	Please enquire.