

Product Datasheet

Rat Transforming Growth Factor Beta 3 (TGFb3) ELISA Kit (orb1208215)

Description	Rat Transforming Growth Factor Beta 3 (TGFb3) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Uniprot ID	Q07258
Sensitivity	5.8pg/mL
Entrez	25717
Expiration Date	Please enquire.