

Product Datasheet

Rat Interferon Regulatory Factor 4 (IRF4) ELISA Kit (orb1208220)

Description	Rat Interferon Regulatory Factor 4 (IRF4) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	D3ZEX6
Sensitivity	0.133ng/mL
Entrez	291100
Expiration Date	Please enquire.