

Product Datasheet

Rat Interleukin 2 Receptor Alpha (IL2Ra) ELISA Kit (orb1208225)

Description	Rat Interleukin 2 Receptor Alpha (IL2Ra) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P26897
Sensitivity	0.051ng/mL
Entrez	25704
Expiration Date	Please enquire.