

Product Datasheet

Rat Interleukin 25 (IL25) ELISA Kit (orb1208226)

Description	Rat Interleukin 25 (IL25) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	D3ZLB1
Sensitivity	32pg/mL
Entrez	501996
Expiration Date	Please enquire.