

Product Datasheet

Rat Insulin Like Growth Factor 1 Receptor (IGF1R) ELISA Kit (orb1208232)

Description	Rat Insulin Like Growth Factor 1 Receptor (IGF1R) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P24062
Sensitivity	0.062ng/mL
Entrez	25718
Expiration Date	Please enquire.