

Product Datasheet

Rat Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit (orb1208239)

Description	Rat Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q5EGZ1
Sensitivity	0.25ng/mL
Entrez	302668
Expiration Date	Please enquire.