

Product Datasheet

Guinea pig Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit (orb1208245)

Description	Guinea pig Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit
Reactivity	Guinea pig
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Z1Y6
Sensitivity	7.6pg/mL
Entrez	100135628
Expiration Date	Please enquire.