

Product Datasheet

Monkey Interleukin 8 (IL8) ELISA Kit (orb1208253)

Description	Monkey Interleukin 8 (IL8) ELISA Kit
Reactivity	Monkey
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P67813
Sensitivity	5.3pg/mL
Entrez	613028
Expiration Date	Please enquire.