

Product Datasheet

Mouse Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit (orb1208254)

Description	Mouse Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P35343
Sensitivity	12.4pg/mL
Entrez	12765
Expiration Date	Please enquire.