

Product Datasheet

Canine Interleukin 2 (IL2) ELISA Kit (orb1208262)

Description	Canine Interleukin 2 (IL2) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q29416
Sensitivity	6.2pg/mL
Entrez	403989
Expiration Date	Please enquire.