

Product Datasheet

Canine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit (orb1208268)

Description	Canine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit
Reactivity	Canine
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9BEH0
Sensitivity	14.35pg/mL
Entrez	403660
Expiration Date	Please enquire.