

Product Datasheet

Guinea pig Interleukin 1 Alpha (IL1a) ELISA Kit (orb1208269)

Description	Guinea pig Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	HOW0J5
Sensitivity	6.6pg/mL
Entrez	100135555
Expiration Date	Please enquire.