

Product Datasheet

Canine Interleukin 18 (IL18) ELISA Kit (orb1208271)

Description	Canine Interleukin 18 (IL18) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XSR0
Sensitivity	5.9pg/mL
Entrez	403796
Expiration Date	Please enquire.