

Product Datasheet

Mouse Interleukin 18 Binding Protein (IL18BP) ELISA Kit (orb1208272)

Description	Mouse Interleukin 18 Binding Protein (IL18BP) ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9Z0M9
Sensitivity	45pg/mL
Entrez	16068
Expiration Date	Please enquire.