

Product Datasheet

Mouse Interleukin 11 Receptor Alpha (IL11Ra) ELISA Kit (orb1208276)

Description	Mouse Interleukin 11 Receptor Alpha (IL11Ra) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q64385
Sensitivity	0.061ng/mL
Entrez	16157
Expiration Date	Please enquire.