

Product Datasheet

Guinea pig Interleukin 10 (IL10) ELISA Kit (orb1208277)

Description	Guinea pig Interleukin 10 (IL10) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	HOVDI2
Sensitivity	5.2pg/mL
Entrez	100715618
Expiration Date	Please enquire.