

Product Datasheet

Mouse Growth Factor, Augmenter Of Liver Regeneration (GFER) ELISA Kit (orb1208287)

Description	Mouse Growth Factor, Augmenter Of Liver Regeneration (GFER) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P56213
Sensitivity	29pg/mL
Entrez	11692
Expiration Date	Please enquire.