

Product Datasheet

Guinea pig Endothelin 1 (EDN1) ELISA Kit (orb1208295)

Description	Guinea pig Endothelin 1 (EDN1) ELISA Kit
Reactivity	Guinea pig
Range	7.81-500pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P97740
Sensitivity	2.82pg/mL
Entrez	100726197
Expiration Date	Please enquire.