

Product Datasheet

Guinea pig Brain Derived Neurotrophic Factor (BDNF) ELISA Kit (orb1208298)

Description	Guinea pig Brain Derived Neurotrophic Factor (BDNF) ELISA Kit
Reactivity	Guinea pig
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	070183
Sensitivity	0.059ng/mL
Entrez	100730257
Expiration Date	Please enquire.