

Product Datasheet

Human Tumor Necrosis Factor Alpha Induced Protein 2 (TNFaIP2) ELISA Kit (orb1208366)

Description	Human Tumor Necrosis Factor Alpha Induced Protein 2 (TNFaIP2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q03169
Sensitivity	0.109ng/mL
Entrez	7127
Expiration Date	Please enquire.