

Product Datasheet

Human Platelet Derived Growth Factor Receptor Like Protein (PDGFRL) ELISA Kit (orb1208370)

Description	Human Platelet Derived Growth Factor Receptor Like Protein (PDGFRL) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q15198
Sensitivity	0.041ng/mL
Entrez	5157
Expiration Date	Please enquire.