

Product Datasheet

Human Platelet Derived Growth Factor Receptor Beta (PDGFRb) ELISA Kit (orb1208371)

Description	Human Platelet Derived Growth Factor Receptor Beta (PDGFRb) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P09619
Sensitivity	0.063ng/mL
Entrez	5159
Expiration Date	Please enquire.