

Product Datasheet

Human Lymphocyte Antigen 96 (LY96) ELISA Kit (orb1208376)

Description	Human Lymphocyte Antigen 96 (LY96) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Y6Y9
Sensitivity	0.060ng/mL
Entrez	23643
Expiration Date	Please enquire.