

Product Datasheet

Human Interleukin 28B (IL28B) ELISA Kit (orb1208385)

Description	Human Interleukin 28B (IL28B) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8IZI9
Sensitivity	5.6pg/mL
Entrez	282617
Expiration Date	Please enquire.