

Product Datasheet

Human Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit (orb1208388)

Description	Human Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q969J5
Sensitivity	0.127ng/mL
Entrez	116379
Expiration Date	Please enquire.