

Product Datasheet

Human Interleukin 1 Epsilon (IL1e) ELISA Kit (orb1208394)

Description	Human Interleukin 1 Epsilon (IL1e) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9UHA7
Sensitivity	5.1pg/mL
Entrez	27179
Expiration Date	Please enquire.