

Product Datasheet

Human Interleukin 11 Receptor Alpha (IL11Ra) ELISA Kit (orb1208398)

Description	Human Interleukin 11 Receptor Alpha (IL11Ra) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q14626
Sensitivity	0.112ng/mL
Entrez	3590
Expiration Date	Please enquire.