

Product Datasheet

Human B-Lymphoid Tyrosine Kinase (BLK) ELISA Kit (orb1208420)

Description	Human B-Lymphoid Tyrosine Kinase (BLK) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51451
Sensitivity	0.071ng/mL
Entrez	640
Expiration Date	Please enquire.