

## Product Datasheet

# Human ATPase, H<sup>+</sup>/K<sup>+</sup> Exchanging Alpha Polypeptide (ATP4a) ELISA Kit (orb1208422)

|                          |                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human ATPase, H <sup>+</sup> /K <sup>+</sup> Exchanging Alpha Polypeptide (ATP4a) ELISA Kit |
| <b>Reactivity</b>        | Human                                                                                       |
| <b>Range</b>             | 0.156-10ng/mL                                                                               |
| <b>Note</b>              | For research use only                                                                       |
| <b>Application notes</b> | sandwich                                                                                    |
| <b>Sample Types</b>      | tissue homogenates, cell lysates or other biological fluids.                                |
| <b>Uniprot ID</b>        | <b>P20648</b>                                                                               |
| <b>Sensitivity</b>       | 0.055ng/mL                                                                                  |
| <b>Entrez</b>            | <b>495</b>                                                                                  |
| <b>Expiration Date</b>   | Please enquire.                                                                             |