

Product Datasheet

Goat Interleukin 6 (IL6) ELISA Kit (orb1208424)

Description	Goat Interleukin 6 (IL6) ELISA Kit
Reactivity	Goat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	E3UT57
Sensitivity	6.1pg/mL
Entrez	100860785
Expiration Date	Please enquire.