

Product Datasheet

Rabbit Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit (orb1208428)

Description	Rabbit Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	5.9pg/mL
Entrez	100008899
Expiration Date	Please enquire.