

Product Datasheet

Rat Transforming Growth Factor Alpha (TGFa) ELISA Kit (orb1208442)

Description	Rat Transforming Growth Factor Alpha (TGFa) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01134
Sensitivity	6.7pg/mL
Entrez	24827
Expiration Date	Please enquire.