

Product Datasheet

Mouse Interferon Regulatory Factor 5 (IRF5) ELISA Kit (orb1208452)

Description	Mouse Interferon Regulatory Factor 5 (IRF5) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P56477
Sensitivity	0.048ng/mL
Entrez	27056
Expiration Date	Please enquire.