

Product Datasheet

Mouse Interferon Regulatory Factor 3 (IRF3) ELISA Kit (orb1208453)

Description	Mouse Interferon Regulatory Factor 3 (IRF3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P70671
Sensitivity	0.048ng/mL
Entrez	54131
Expiration Date	Please enquire.