

Product Datasheet

Rat Interleukin 7 (IL7) ELISA Kit (orb1208460)

Description	Rat Interleukin 7 (IL7) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56478
Sensitivity	6.5pg/mL
Entrez	25647
Expiration Date	Please enquire.