

Product Datasheet

Rabbit Interleukin 2 (IL2) ELISA Kit (orb1208472)

Description	Rabbit Interleukin 2 (IL2) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	077620
Sensitivity	5.7pg/mL
Entrez	100302458
Expiration Date	Please enquire.