

Product Datasheet

Rabbit Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit (orb1208480)

Description	Rabbit Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit
Reactivity	Rabbit
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P26890
Sensitivity	12.8pg/mL
Entrez	100009558
Expiration Date	Please enquire.