

Product Datasheet

Chicken Interleukin 1 Beta (IL1b) ELISA Kit (orb1208482)

Description	Chicken Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Gallus
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	073909
Sensitivity	5.7pg/mL
Entrez	395196
Expiration Date	Please enquire.