

Product Datasheet

Mouse Interleukin 19 (IL19) ELISA Kit (orb1208485)

Description	Mouse Interleukin 19 (IL19) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8CJ70
Sensitivity	6.3pg/mL
Entrez	329244
Expiration Date	Please enquire.