

Product Datasheet

Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit (orb1208494)

Description	Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O88786
Sensitivity	0.126ng/mL
Entrez	16165
Expiration Date	Please enquire.