

Product Datasheet

Chicken Interleukin 13 (IL13) ELISA Kit (orb1208496)

Description	Chicken Interleukin 13 (IL13) ELISA Kit
Reactivity	Gallus
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5W4U0
Sensitivity	5.3pg/mL
Entrez	474361
Expiration Date	Please enquire.