

Product Datasheet

Mouse Insulin Like Growth Factor Binding Protein 7 (IGFBP7) ELISA Kit (orb1208505)

Description	Mouse Insulin Like Growth Factor Binding Protein 7 (IGFBP7) ELISA Kit
Reactivity	Mouse
Range	6.25-400pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q61581
Sensitivity	2.34pg/mL
Entrez	29817
Expiration Date	Please enquire.